

Cervine

Anterior Cervical Fusion Solution

- ◎ One-step Locking
- ◎ 2.0mm Low Profile
- ◎ Wide Graft Visualization



SIMPLICITY
FLEXIBILITY
RELIABILITY




Anterior Cervical Spine System

Diverse Screw Options

Provides intraoperative flexibility to better fit the patient's anatomy



Fixed Screw



Variable Screw



Self-drilling



Self-tapping

Optimal Implant Design

- One-level, two-level, three-level plates available, length from 23mm – 76mm
- Low profile minimizes the soft tissue irritation



One-level Plate Two-level Plate Three-level plate

More Efficient Anterior Cervical Spine System




Midline Alignment Notch

Visual lock verification

Graft visualization

Fixed screw/
Variable screw



Low profile

2mm

