

**KANOHUI**

A  **Medtronic** Company



# **MonoLoc**

**Locking Plates 1.5/2.0mm**

### Indication

#### Possible indications for implants of size 1.5 include:

- Fractures, osteotomies and arthrodeses of the phalanges and metacarpals
- Replantations and reconstructions of the phalanges and metacarpals

#### Possible indications for implants of size 2.0 include:

- Fractures of the phalanges, metacarpals, wrist bones and distal radius (double-plate technique)
- Osteotomies and arthrodeses of the interphalangeal joints

### Features & Benefits

- Coaxial locking hole accept both locking screw and cortex screw in the same round conical plate hole;




- Low profile plates and fully countersunk screws to minimize soft tissue irritation;



- Anatomically precontoured plate minimizes need for bending;



- A fixed-angle construct is particularly advantageous in osteoporotic bone or in fractures near the joint;



- Self-tapping screws with low profile head allows easy insertion and reduce irritation to the soft tissue;



CHANGZHOU KANGHUI MEDICAL INNOVATION CO., LTD

Add: No.11, North Changjiang Road, Xinbei Zone, Changzhou, Jiangsu 213022, P.R.China  
Tel: +86-519-85139851 Fax: +86-519-85128628  
www.kanghui.com  
Shanghai Office  
Add: Building 18, No.1000, Jinhai Road, PuDong District, Shanghai, China 201206  
Tel: +86-21-50319916 Fax: +86-21-50312913